
**KABO
TECH®**

Waterbaden

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All rights reserved. No part of this publication may be copied or published by means of printing, photocopying, microfilm or otherwise without prior written consent of LaboTech®. This restriction also applies to the corresponding drawings and diagrams.

LaboTech® reserves the right to change parts of the machine at any time without prior or direct notice to the client. The contents of this publication are subject to change without notice.

This publication is to be used for the standard version of the machine only. Thus LaboTech® cannot be held responsible for any damage resulting from the application of this publication to the version actually delivered to you.

This publication has been written with great care. However, LaboTech® cannot be held responsible for any errors occurring in this publication or for their consequences.

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TABLE OF CONTENTS

Technical data	19
Safety instructions	20
Package contents	21
Description	21
Electrical safety	21
Mains plug replacement (U.K. & Ireland only)	22
Using an extension cable	22
Instructions for use	22
Maintenance	23
Possible errors	24
LaboTech® after-sales service	24
Guarantee	24
Unwanted apparatus and the environment	24

The following symbols are used throughout this manual:

Denotes risk of personal injury, loss of life or damage to the apparatus in case of nonobservance of the instructions in this manual.

Avoid bodily contact with the waterbath.

Fire hazard.
Keep inflammable products away from the waterbath.

Ambient conditions.

For indoor use and for laboratory applications only.

Do not use in explosive environments.

Geographical height

(max) m 2.000

Ambient temperatures °C 5 - 40

Max. r.h. from 80 % at 31 °C and decreasing to 50%at40°C.

Technical data	OWB14N	DWB14N (180 mm)	DWB16N	DWB16N (180 mm)	DWB23N	DWB23N (180mm)
Voltage (V)	220 ±10%	220 ±10%	220 ±10%	220 ±10%	220 ±10%	220 ±10%
Frequency (Hz)	50	50	50	50	50	50
Current (A)	10	10	10	10	10	10
Power (W)	2200	2200	2200	2200	2200	2200
Heatingtemp. (max °C)	120	120	120	120	120	120
Dimension (lxbxh) (mm)	605x290x 195	605x290x 285	905x290x 195	905x290x 285	455x435x 195	455x435x 285
Bathdimension (lxbxh) (mm)	600x170x 100	600x170x 180	900x170x 100	900x170x 180	450x320x 100	450x320x 180
Weight (Kg)	10,4	12,0	15,3	17,4	11,8	13,4
Ordernumber	11.14.0001	11.14.0181	11.16.0001	11.16.0181	11.23.0001	11.23.0181

Safety instructions.

When using electrical apparatus, always observe the safety regulations applicable in your country to reduce the risk of fire, electric shock and personal injury. Read the following safety instructions before attempting to operate this product.

Keep these instructions in a safe place.

General..

1 Keep work area clean.

Cluttered areas and benches can cause accidents.

2 Consider work area environment.

Do not expose electrical apparatus to humidity. Keep work area well lit. Do not use electrical apparatus in the presence of inflammable liquids or gases.

3 Guard against electric shock.

Prevent body contact with earthed surfaces (e.g. pipes, radiators, cookers and refrigerators).

For use under extreme conditions (e.g. high humidity) electric safety can be improved by inserting an isolating transformer or a (FI) earth-leakage circuit-breaker.

4 Keep unauthorized persons away.

Do not let unauthorized persons touch the apparatus or the mains cable. Keep unauthorized people away from the work area.

5 Dress properly.

Do not wear loose clothing or jewellery. Wear protective hair covering to keep long hair out of the way.

6 Wear safety goggles.

Always wear safety goggles in view of possible spattering of the product.

7 Do not overreach.

Keep proper footing and balance at all times.

8 Stay alert.

Watch what you are doing. Use common sense. Do not operate the apparatus when you are tired.

9 Do not abuse cord.

Never pull the cord to disconnect from the socket. Keep the cord away from heat, oil and sharp edges.

10 Maintain the apparatus with care.

Keep the apparatus in good condition and clean for better and safer performance. Follow the instructions for maintenance and changing accessories. Inspect the apparatus cords at regular intervals and if damaged, have them repaired by a LaboTech® authorized repair agent.

Inspect the extension cords periodically and replace them if damaged. Keep all controls dry, clean and free from oil and grease.

11 Check for damaged parts.

Before using the apparatus, carefully check it for damage to ensure that it will operate properly and perform its intended function. Have damaged guards or other defective parts repaired or replaced as instructed. Do not use the apparatus if the switch is defective. Have the switch replaced by a LaboTech® authorized repair agent.

12 Have your apparatus repaired by a LaboTech® authorized repair agent.

This electrical apparatus is in accordance with the relevant safety regulations. To avoid danger, electric appliance must only be repaired by qualified technicians.




fig.A

fig.c

fig.B

Package contents.

The package contains:

- 1 waterbath
- 1 cover
- 1 instruction manual
- 1 Guarantee card

Optional:

- External temperature sensor
- Ringset
- Constant level
- Standassemblyset
- Check for damage to the apparatus, parts or accessoires which may have occurred during transport.
- Take the time to thoroughly read and understand this manual prior to operation.

Description (fig. A, B and C).

Your LaboTech® waterbath has been designed for professional use. The heat is evenly distributed over the entire surface.

Thanks to the newly developed temperature regulation, the heat can be set very accurately.

- 1 ON/OFF-switch with indicator
- 2 Connector for external temperature sensor
- 3 Cistern
- 4 Display
- 5 Touch key with LED ACT. TEMP.
- 6 Touch key with LED SET TEMP.
- 7 Touch key with LED 100 % HEATING
- 8 Touch key with LED EXT. SENSOR
- 9 Touch key with LED INT. SENSOR
- 10 Touch key with LED 50 % HEATING
- 11 LED SET TEMP.
- 12 Touch key with LED 25 % HEATING
- 13 Touch key with LED STAND BY
- 14 Reset button

Electrical safety.

The apparatus has been designed for one voltage only. Always check that the power supply corresponds to the voltage on the rating plate.


fig.D

Mains plug replacement (U.K. & Ireland only) (Fig. D)

- Should your mains plug need replacing and you are competent to do this, proceed as instructed below. If you are in doubt, contact a LaboTech® authorized repair agent or a qualified electrician.
- Disconnect the plug from the supply.
- Cut off the plug and dispose of it safely; a plug with bared copper conductors is dangerous if engaged in a live socket outlet.
- Only fit 13 Amperes BS1363A approved plugs fitted with the correctly rated fuse(1).
- The cable wire colours, or a letter, will be marked at the connection points of most good quality plugs. Attach the wires to their respective points in the plug (see below). Brown is for Live (L) (2) and Blue is for Neutral (N) (4).
- Before replacing the top cover of the mains plug ensure that the cable restraint(3) is holding the outer sheath of the cable firmly and that the two leads are correctly fixed at the terminal screws.

Never use a light socket.
Never connect the live (L) or neutral (N) wires to the earth pin marked E or.

Using an extension cable.

If an extension cable is required, use an approved triple-core extension cable suitable for the power input of this apparatus (see technical data). The minimum conductor size is 2.5 mm².

Instructions for use.

Always observe the safety instructions and applicable regulations.

Take care of the safety instructions.

Never use the apparatus to heat lightly inflammable liquids.

Wear isolating gloves and/or use pilers to manipulate the heated products.

Do not empty the waterbath in heated condition.

- Use grounded wall sockets only.
- Fill up the waterbath.

Fill up (in liters)	DWB 14N min./max.	DWB16N min./max.	DWB23N min./max.
	3,5/8,0	5,3/11,9	5,0/11,2
180 mm	3,5/17,0	5,3/25,7	5,0/24,19

- Set the ON/OFF-switch (1) to ON-position.
The waterbath will go to 'stand by'-mode.

Using an internal temperature sensor,

- Press the touch key INT. SENSOR(9). The LED lights.

Option:

Using the external temperature sensor.

- Set the ON/OFF-switch(1) to OFF-position.
- Insert the plug of the external temperature sensor into the connector(2).
- Place the temperature sensor in the product to be heated.
- Set the ON/OFF-switch(1) to ON-position.
- Press the touch key EXT. SENSOR(8).
The LED lights.
The waterbath may be hotter than the temperature indicated by the external temperature sensor (heat transfer losses).
For the correct temperature, please refer to the section "Using the internal temperature sensor".

Trouble shooting:

Display shows 458.

Cause: Temperature sensor failure

Waterbath doesn't heat up.

Cause: Element protection active

Check waterlevel (see instructions for use).

Push reset button (see fig. C) until you hear a "click". (First remove black nut).

Setting the required temperature.

- Press the touch key SET TEMP.(6). The LED lights.
- Rotate the SET TEMP.(11) wheel to set the required temperature.

The display indicates the set temperature.

Adjusting the heating power.

The value of the heating power depends on the size of the product to be heated and/or the set temperature.

Setting to 25 % of the heating power:

- Press the touch key 25 % HEATING(12).
The LED lights.

Setting to 50 % of the heating power:

- Press the touch key 50 % HEATING(10).
The LED lights.

Setting to full power:

- Press the touch key 100 % HEATING(7).
The LED lights.

Interruptions.

- Press the touch key STAND BY(13). The LED lights.

Switching off.

- Set the ON/OFF switch(1) to OFF-position.

The waterbath will cool only very slowly.

Maintenance.

Your LaboTech® Apparatus has been designed to operate over a long period of time with a minimum of maintenance. Continuous satisfactory operation depends upon proper care and regular cleaning.

Cleaning.

Regularly clean the housing with a soft cloth.

The waterbath can be very hot.

Labotech After Sales service.

All LaboTech® apparatus are thoroughly tested before leaving the factory. However, if the apparatus needs repair, please contact your dealer or your LaboTech® authorized repair agent.

Guarantee.

This apparatus is guaranteed by the manufacturer in accordance with the terms of the enclosed guarantee card.

Unwanted apparatus and the environment.

Take your apparatus to a LaboTech® authorized repair agent, where it will be disposed of in an environmentally safe way.